



**PLEX Elite 9000 from ZEISS**  
Swept-Source OCT





# Uncovering the undiscovered.

ZEISS PLEX Elite 9000

// INNOVATION  
MADE BY ZEISS


Ultra-wide angiography *En face* montage  
Image courtesy of Prof. G. Querques, San Raffaele Ospedale, Milan, Italy

A new idea is often the start of scientific discovery. But it is transformational new technology that often enables researchers to act upon these ideas and to explore previously unreachable frontiers. PLEX® Elite 9000 from ZEISS is just such a technology. By inviting researchers into a new world of structural and microvascular clarity of the anatomy, PLEX Elite 9000 is foundational to the future of retina research and to the understanding of the development of retinal disease.

## ZEISS PLEX Elite 9000

- **SEE** deeper, wider and in more detail
- **STUDY** early mechanisms of micro- and neovascularization of the posterior segment from vitreous to sclera
- **EXPLORE** the progression of retinal and choroidal pathology, such as CNV
- **IMPROVE** understanding of choroid physiopathology
- **EVALUATE** the mechanism of retina and choroid response to a therapy






HD Spotlight 16 mm B-scan of choroidal excavation  
Image courtesy of Prof. F.G .Holz, Universitäts Augenklinik, Bonn, Germany

## Explore deeper meanings


ZEISS PLEX Elite 9000 Swept-Source OCT allows clinical researchers the potential to see deeper, wider and in more detail from the vitreous to the sclera in the posterior segment.

**NEW Montage** scan acquisition workflow and export features showcase the ability to rapidly acquire an ultra-wide OCT angiography *En face* montage for unprecedented visualization of retinal vasculature with a field of view up to 70°





**NEW Ultra-wide 15x9** high-density scan reveals the widest field of view captured in a single OCT angiography scan




Ultra-widefield AngioPlex montage of proliferative diabetic retinopathy, superficial layer  
Image courtesy of Prof Korobelnik, University Hospital Pellegrin, Bordeaux, France



Ultra-wide 15x9 AngioPlex map of choroidal neovascular membrane, RPE/RPE Fit layer  
Image courtesy of Prof. Rosenfeld, MD, Bascom Palmer Eye Institute, Miami, FL



On-the-fly B-scan through the atrophic area from the 12x12 OCT En face cube



12x12 AngioPlex map, superficial layer and corresponding B-scan of a vitreomacular traction  
 Images courtesy of Prof. Rosenfeld, MD, Bascom Palmer Eye Institute, Miami, FL

**Standard view** AngioPlex® maps provide a full view of the retina at 3x3, 6x6, 9x9 or 12x12

**HD spotlight** high-detail B-scan up to 16 mm

**On-the-fly** B-scan to define a custom slice through any angle

**UHD (ultra-high definition)** cubes with excellent image quality of fundus, clear visualization of vasculature and ocular structures at any depth from vitreous to sclera


**AngioPlex OCT** angiography for ultra-clear 3D microvascular visualizations powered by OMAG<sup>c</sup>

**FastTrac™** real-time tracking of eye movement for motion-artifact compensation

**Up-to-date software and hardware technology over two years** to keep you at the cutting edge of development in the area of Swept-Source OCT

# A R I Network

Discover. Collaborate. Understand.



The Advanced Retina Imaging (A R I) Network, with ZEISS PLEX Elite 9000 at its core, brings together the **expertise of leading clinicians and researchers** around the world with scientists and developers **at ZEISS** to accelerate the development of innovations to benefit patients today and in the future.

Through an active exchange of ideas and findings, the aim of the A R I Network is to drive the development of new clinical applications and future OCT technologies.

[www.zeiss.com/arinetwork](http://www.zeiss.com/arinetwork)

# Technical Specifications

## OCT Imaging

Methodology	Swept-Source OCT
Optical source	Swept-Source tunable laser: center wavelength between 1040 nm and 1060 nm
Scan speed	100,000 A-scans/sec
A-scan depth	3.0 mm (in tissue)
Axial resolution (optical)	6.3 $\mu\text{m}$ (in tissue)
Axial resolution (digital)	1.95 $\mu\text{m}$ (in tissue)
Transverse resolution*	20 $\mu\text{m}$ (*transverse [Lateral] resolution is calculated from the beam size at the pupil)
Field of view	56°
Minimum pupil diameter	2.5 mm

## Fundus imaging

Methodology	Line-scanning ophthalmoscope (LSO) – live fundus image during alignment and during OCT scan
Optical source	Super-luminescent diode (SLD) 750 nm
Field of view	36° W x 30° H
Frame rate	> 20 Hz

## Iris imaging

Methodology	CCD camera
Resolution	1280x1024

ZEISS PLEX Elite 9000 is CE Marked and is available for sale in selected countries and in the United States.

**CE** 0297



**Carl Zeiss Meditec, Inc.**  
5160 Hacienda Drive  
Dublin, CA 94568  
USA  
[www.zeiss.com/med](http://www.zeiss.com/med)  
[www.zeiss.com/med/contacts](http://www.zeiss.com/med/contacts)



**Carl Zeiss Meditec AG**  
Goeschwitzer Str. 51-52  
07745 Jena  
Germany  
[www.zeiss.com/med](http://www.zeiss.com/med)  
[www.zeiss.com/med/contacts](http://www.zeiss.com/med/contacts)

**EN\_31\_020\_0001V / US\_31\_020\_0001V** Printed in the United States. CZ-IX/2017 International edition: Only for sale in selected countries.  
The contents of the brochure may differ from the current status of approval of the product or service offering in your country. Please contact our regional representative for more information. Subject to changes in design and scope of delivery and as a result of ongoing technical development. PLEX Elite, AngioPlex and FastTrac are either trademarks or registered trademarks of Carl Zeiss Meditec AG or other companies of the ZEISS Group in Germany and/or other countries.  
© Carl Zeiss Meditec, Inc. 2017. All rights reserved.